

Anti-Cd274/PD-L1 antibody [29E.2A3] (STJA0001825)

STJA0001825

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-Programmed cell death 1 ligand 1 is suitable for use in Flow Cytometry and Immunohistochemistry research applications.
Applications	FC, IHC
Host/Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	29E.2A3
Concentration	0.5 mg/mL
Conjugation	Unconjugated
Purification	
Dilution Range	
Formulation	PBS with 0.05% Proclin300
Isotype	IgG2bk
Storage	Recommend storing at 2-8°C. Do not freeze this product.
Instruction	

TARGET INFORMATION

Gene ID	29126
Gene Symbol	CD274
Uniprot ID	PD1L1_HUMAN
Immunogen	
Immunogen	
Region	
Specificity	
Immunogen	
Sequence	