

Anti-CD274/PD-L1 antibody [29E.2A3] {PerCP/Cyanine5.5} (STJA0001829)

STJA0001829

GENERAL INFORMATION

| | |
|--------------------------|---|
| Product Type | Primary antibodies |
| Short Description | Mouse monoclonal antibody anti-CD274/PD-L1 is suitable for use in Flow Cytometry research applications. |
| Applications | FC |
| Host/Source | Mouse |
| Reactivity | Human |

PRODUCT PROPERTIES

| | |
|----------------------------|--|
| Clonality | Monoclonal |
| Clone ID | 29E.2A3 |
| Concentration | |
| Conjugation | PerCP/Cyanine5.5 |
| Purification | |
| Dilution Range | |
| Formulation | PBS with 0.05% Proclin300, 1% BSA |
| Isotype | IgG2bk |
| Storage Instruction | Recommend storing between 2-8°C and protecting from prolonged exposure to light. Do not freeze this product. |

TARGET INFORMATION

| | |
|---------------------------|-----------------------------|
| Gene ID | 29126 |
| Gene Symbol | CD274 |
| Uniprot ID | PD1L1_HUMAN |
| Immunogen | |
| Immunogen Region | |
| Specificity | |
| Immunogen Sequence | |