

Anti-CD196/CCR6 antibody [G034E3] {PE/Cyanine7} (STJA0001915)

STJA0001915

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-C-C chemokine receptor type 6 is suitable for use in Flow Cytometry research applications.
Applications	FC
Host/Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	G034E3
Concentration	
Conjugation	PE/Cyanine7
Purification	
Dilution Range	
Formulation	PBS with 0.05% Proclin300, 1% BSA
Isotype	IgG2bk
Storage Instruction	Recommend storing between 2-8°C and protecting from prolonged exposure to light. Do not freeze this product.

TARGET INFORMATION

Gene ID	1235
Gene Symbol	CCR6
Uniprot ID	CCR6_HUMAN
Immunogen	
Immunogen Region	
Specificity	
Immunogen Sequence	