

Anti-CD197/CCR7 antibody [G043H7] {PE/Cyanine5} (STJA0001916)

STJA0001916

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-C-C chemokine receptor type 7 is suitable for use in Flow Cytometry research applications.
Applications	FC
Host/Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	G043H7
Concentration	
Conjugation	PE/Cyanine5
Purification	
Dilution Range	
Formulation	PBS with 0.05% Proclin300, 1% BSA
Isotype	IgG2ak
Storage Instruction	Recommend storing between 2-8°C and protecting from prolonged exposure to light. Do not freeze this product.

TARGET INFORMATION

Gene ID	1236
Gene Symbol	CCR7
Uniprot ID	CCR7_HUMAN
Immunogen	
Immunogen Region	
Specificity	
Immunogen Sequence	