

Anti-CD200 antibody [OX-104] {APC} (STJA0001921)

STJA0001921

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-OX-2 membrane glycoprotein is suitable for use in Flow Cytometry research applications.
Applications	FC
Host/Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	OX-104
Concentration	
Conjugation	APC
Purification	
Dilution Range	
Formulation	PBS with 0.05% Proclin300, 1% BSA
Isotype	IgG1k
Storage Instruction	Recommend storing between 2-8°C and protecting from prolonged exposure to light. Do not freeze this product.

TARGET INFORMATION

Gene ID	4345
Gene Symbol	CD200
Uniprot ID	OX2G_HUMAN
Immunogen	
Immunogen Region	
Specificity	
Immunogen Sequence	