

## Anti-CD206-antibody [15.2] {FITC} (STJA0001922)

STJA0001922

### GENERAL INFORMATION

<b>Product Type</b>	Primary antibodies
<b>Short Description</b>	Mouse monoclonal antibody anti-CD206 is suitable for use in Flow Cytometry.
<b>Applications</b>	FC
<b>Host/Source</b>	Mouse
<b>Reactivity</b>	Human

### PRODUCT PROPERTIES

<b>Clonality</b>	Monoclonal
<b>Clone ID</b>	15.2
<b>Concentration</b>	
<b>Conjugation</b>	FITC
<b>Purification</b>	
<b>Dilution Range</b>	
<b>Formulation</b>	PBS with 0.05% Proclin300, 1% BSA
<b>Isotype</b>	IgG1, k
<b>Storage Instruction</b>	Store at 2-8°C and protected from prolonged exposure to light. Do not freeze.

### TARGET INFORMATION

<b>Gene ID</b>	4360
<b>Gene Symbol</b>	MRC1
<b>Uniprot ID</b>	MRC1_HUMAN
<b>Immunogen</b>	
<b>Immunogen Region</b>	
<b>Specificity</b>	
<b>Immunogen Sequence</b>	