

Anti-CD39/ENTPD1 antibody [A1] {AF488} (STJA0001934)

STJA0001934

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-Ectonucleoside triphosphate diphosphohydrolase 1 is suitable for use in Flow Cytometry research applications.
Applications	FC
Host/Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	A1
Concentration	
Conjugation	AF488
Purification	
Dilution Range	
Formulation	PBS with 0.05% Proclin300, 1% BSA
Isotype	IgG1k
Storage Instruction	Recommend storing between 2-8°C and protecting from prolonged exposure to light. Do not freeze this product.

TARGET INFORMATION

Gene ID	953
Gene Symbol	ENTPD1
Uniprot ID	ENTP1_HUMAN
Immunogen	
Immunogen	
Region	
Specificity	
Immunogen	
Sequence	