

## Anti-CD61/ITGB3 antibody [VI-PL2] {APC} (STJA0001938)

STJA0001938

### GENERAL INFORMATION

|                          |   |
|--------------------------|---|
| <b>Product Type</b>      | Primary antibodies  |
| <b>Short Description</b> | Mouse monoclonal antibody anti-Integrin beta-3 is suitable for use in Flow Cytometry research applications. |
| <b>Applications</b>      | FC  |
| <b>Host/Source</b>       | Mouse   |
| <b>Reactivity</b>        | Human   |

### PRODUCT PROPERTIES

|                            |  |
|----------------------------|--|
| <b>Clonality</b>           | Monoclonal   |
| <b>Clone ID</b>            | VI-PL2   |
| <b>Concentration</b>       |  |
| <b>Conjugation</b>         | APC  |
| <b>Purification</b>        |  |
| <b>Dilution Range</b>      |  |
| <b>Formulation</b>         | PBS with 0.05% Proclin300, 1% BSA  |
| <b>Isotype</b>             | IgG1k  |
| <b>Storage Instruction</b> | Recommend storing between 2-8°C and protecting from prolonged exposure to light. Do not freeze this product. |

### TARGET INFORMATION

|                           |                            |
|---------------------------|----------------------------|
| <b>Gene ID</b>            | <a href="#">3690</a>       |
| <b>Gene Symbol</b>        | <a href="#">ITGB3</a>      |
| <b>Uniprot ID</b>         | <a href="#">ITB3_HUMAN</a> |
| <b>Immunogen</b>          |                            |
| <b>Immunogen Region</b>   |                            |
| <b>Specificity</b>        |                            |
| <b>Immunogen Sequence</b> |                            |