

Anti-Cd14 antibody [Sa14-2] {FITC} (STJA0001971)

STJA0001971

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rat monoclonal antibody anti-Monocyte differentiation antigen CD14 is suitable for use in Flow Cytometry research applications.
Applications	FC
Host/Source	Rat
Reactivity	Mouse

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	Sa14-2
Concentration	0.5 mg/mL
Conjugation	FITC
Purification	
Dilution Range	
Formulation	PBS with 0.05% Proclin300, 1% BSA
Isotype	IgG2ak
Storage Instruction	Recommend storing between 2-8°C and protecting from prolonged exposure to light. Do not freeze this product.

TARGET INFORMATION

Gene ID	12475
Gene Symbol	Cd14
Uniprot ID	CD14_MOUSE
Immunogen	
Immunogen Region	
Specificity	
Immunogen Sequence	