

Anti-CD23/Fc ϵ 2 antibody [B3B4] {FITC} (STJA0001979)

STJA0001979

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rat monoclonal antibody anti-Low affinity immunoglobulin epsilon Fc receptor is suitable for use in Flow Cytometry research applications.
Applications	FC
Host/Source	Rat
Reactivity	Mouse

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	B3B4
Concentration	
Conjugation	FITC
Purification	
Dilution Range	
Formulation	PBS with 0.05% Proclin300, 1% BSA
Isotype	IgG2ak
Storage Instruction	Recommend storing between 2-8°C and protecting from prolonged exposure to light. Do not freeze this product.

TARGET INFORMATION

Gene ID	14128
Gene Symbol	Fc ϵ 2
Uniprot ID	FCER2_MOUSE
Immunogen	
Immunogen Region	
Specificity	
Immunogen	
Sequence	