

Anti-CD31/Pecam1 antibody [390] {AF647} (STJA0001995)

STJA0001995

GENERAL INFORMATION

| | |
|--------------------------|---|
| Product Type | Primary antibodies |
| Short Description | Rat monoclonal antibody anti-Platelet endothelial cell adhesion molecule is suitable for use in Flow Cytometry research applications. |
| Applications | FC |
| Host/Source | Rat |
| Reactivity | Mouse |

PRODUCT PROPERTIES

| | |
|----------------------------|--|
| Clonality | Monoclonal |
| Clone ID | 390 |
| Concentration | |
| Conjugation | AF647 |
| Purification | |
| Dilution Range | |
| Formulation | PBS with 0.05% Proclin300, 1% BSA |
| Isotype | IgG2ak |
| Storage Instruction | Recommend storing between 2-8°C and protecting from prolonged exposure to light. Do not freeze this product. |

TARGET INFORMATION

| | |
|---------------------------|-----------------------------|
| Gene ID | 18613 |
| Gene Symbol | Pecam1 |
| Uniprot ID | PECA1_MOUSE |
| Immunogen | |
| Immunogen Region | |
| Specificity | |
| Immunogen Sequence | |