

Anti-CD335/Ncr1 antibody [29A1.4] {FITC} (STJA0002010)

STJA0002010

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rat monoclonal antibody anti-Natural cytotoxicity triggering receptor 1 is suitable for use in Flow Cytometry research applications.
Applications	FC
Host/Source	Rat
Reactivity	Mouse

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	29A1.4
Concentration	0.5 mg/mL
Conjugation	FITC
Purification	
Dilution Range	
Formulation	PBS with 0.05% Proclin300, 1% BSA
Isotype	IgG2ak
Storage Instruction	Recommend storing between 2-8°C and protecting from prolonged exposure to light. Do not freeze this product.

TARGET INFORMATION

Gene ID	17086
Gene Symbol	Ncr1
Uniprot ID	NCTR1_MOUSE
Immunogen	
Immunogen Region	
Specificity	
Immunogen Sequence	