

Anti-CD64/Fcgr1 antibody [X54-5/7.1] {FITC} (STJA0002034)

STJA0002034

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-High affinity immunoglobulin gamma Fc receptor I is suitable for use in Flow Cytometry research applications.
Applications	FC
Host/Source	Mouse
Reactivity	Mouse

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	X54-5/7.1
Concentration	0.5 mg/mL
Conjugation	FITC
Purification	
Dilution Range	
Formulation	PBS with 0.05% Proclin300, 1% BSA
Isotype	IgG1k
Storage Instruction	Recommend storing between 2-8°C and protecting from prolonged exposure to light. Do not freeze this product.

TARGET INFORMATION

Gene ID	14129
Gene Symbol	Fcgr1
Uniprot ID	FCGR1_MOUSE
Immunogen	
Immunogen Region	
Specificity	
Immunogen	
Sequence	