

Anti-Cd69 antibody [H1.2F3] {FITC} (STJA0002039)

STJA0002039

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Armenian Hamster monoclonal antibody anti-Early activation antigen CD69 is suitable for use in Flow Cytometry research applications.
Applications	FC
Host/Source	Armenian Hamster
Reactivity	Mouse

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	H1.2F3
Concentration	0.5 mg/mL
Conjugation	FITC
Purification	
Dilution Range	
Formulation	PBS with 0.05% Proclin300, 1% BSA
Isotype	IgG
Storage Instruction	Recommend storing between 2-8°C and protecting from prolonged exposure to light. Do not freeze this product.

TARGET INFORMATION

Gene ID	12515
Gene Symbol	Cd69
Uniprot ID	CD69_MOUSE
Immunogen	
Immunogen Region	
Specificity	
Immunogen Sequence	