

Anti-Fc ϵ 1a antibody [MAR-1] {FITC} (STJA0002045)

STJA0002045

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Armenian Hamster monoclonal antibody anti-High affinity immunoglobulin epsilon receptor subunit alpha is suitable for use in Flow Cytometry research applications.
Applications	FC
Host/Source	Armenian Hamster
Reactivity	Mouse

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	MAR-1
Concentration	0.5 mg/mL
Conjugation	FITC
Purification	
Dilution Range	
Formulation	PBS with 0.05% Proclin300, 1% BSA
Isotype	IgG
Storage	Recommend storing between 2-8°C and protecting from prolonged exposure to light. Do not freeze this product.
Instruction	

TARGET INFORMATION

Gene ID	14125
Gene Symbol	Fc ϵ 1a
Uniprot ID	FCERA_MOUSE
Immunogen	
Immunogen	
Region	
Specificity	
Immunogen	
Sequence	