

Anti-Fc ϵ 1a antibody [MAR-1] {PE} (STJA0002046)

STJA0002046

GENERAL INFORMATION

| | |
|--------------------------|--|
| Product Type | Primary antibodies |
| Short Description | Armenian Hamster monoclonal antibody anti-High affinity immunoglobulin epsilon receptor subunit alpha is suitable for use in Flow Cytometry research applications. |
| Applications | FC |
| Host/Source | Armenian Hamster |
| Reactivity | Mouse |

PRODUCT PROPERTIES

| | |
|----------------------------|--|
| Clonality | Monoclonal |
| Clone ID | MAR-1 |
| Concentration | 0.2 mg/mL |
| Conjugation | PE |
| Purification | |
| Dilution Range | |
| Formulation | PBS with 0.05% Proclin300, 1% BSA |
| Isotype | IgG |
| Storage Instruction | Recommend storing between 2-8°C and protecting from prolonged exposure to light. Do not freeze this product. |

TARGET INFORMATION

| | |
|--------------------|------------------|
| Gene ID | 14125 |
| Gene Symbol | Fc ϵ 1a |
| Uniprot ID | FCERA_MOUSE |
| Immunogen | |
| Region | |
| Specificity | |
| Immunogen | |
| Sequence | |