

Anti-IgD antibody [11-26c.2a] {AF488} (STJA0002060)

STJA0002060

GENERAL INFORMATION

| | |
|--------------------------|---------------------------------------------------------------------------------------------------------------------------------|
| Product Type | Primary antibodies |
| Short Description | Rat monoclonal antibody anti-Ig delta chain C region secreted form is suitable for use in Flow Cytometry research applications. |
| Applications | FC |
| Host/Source | Rat |
| Reactivity | Mouse |

PRODUCT PROPERTIES

| | |
|----------------------------|--------------------------------------------------------------------------------------------------------------|
| Clonality | Monoclonal |
| Clone ID | 11-26c.2a |
| Concentration | 0.5 mg/mL |
| Conjugation | AF488 |
| Purification | |
| Dilution Range | |
| Formulation | PBS with 0.05% Proclin300, 1% BSA |
| Isotype | IgG2ak |
| Storage Instruction | Recommend storing between 2-8°C and protecting from prolonged exposure to light. Do not freeze this product. |

TARGET INFORMATION

| | |
|---------------------------|----------------------------|
| Gene ID | NA |
| Gene Symbol | NA |
| Uniprot ID | IGHD_MOUSE |
| Immunogen | |
| Immunogen Region | |
| Specificity | |
| Immunogen Sequence | |