

Anti-Ly6a antibody [D7] {PE/Cyanine5.5} (STJA0002091)

STJA0002091

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rat monoclonal antibody anti-Lymphocyte antigen 6A-2/6E-1 is suitable for use in Flow Cytometry research applications.
Applications	FC
Host/Source	Rat
Reactivity	Mouse

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	D7
Concentration	0.2 mg/mL
Conjugation	PE/Cyanine5.5
Purification	
Dilution Range	
Formulation	PBS with 0.05% Proclin300, 1% BSA
Isotype	IgG2ak
Storage Instruction	Recommend storing between 2-8°C and protecting from prolonged exposure to light. Do not freeze this product.

TARGET INFORMATION

Gene ID	110454
Gene Symbol	Ly6a
Uniprot ID	LY6A_MOUSE
Immunogen	
Immunogen Region	
Specificity	
Immunogen Sequence	