

Anti-CD366/HAVCR2 antibody [RMT3-23] (STJA0002094)

STJA0002094

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rat monoclonal antibody anti-Hepatitis A virus cellular receptor 2 is suitable for use in Flow Cytometry and Universal research applications.
Applications	FC, Block
Host/Source	Rat
Reactivity	Mouse

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	RMT3-23
Concentration	0.5 mg/mL
Conjugation	Unconjugated
Purification	
Dilution Range	
Formulation	0.2 µm filtered in PBS, pH7.2. Azide Free (AF) /Low Endotoxin (LE) : Contains no stabilizers or preservatives. Endotoxin level is < 2 EU/mg as determined by LAL gel clotting assay.
Isotype	IgG2ak
Storage Instruction	Recommend storing at 2-8°C. Do not freeze this product.

TARGET INFORMATION

Gene ID	84868
Gene Symbol	HAVCR2
Uniprot ID	HAVR2_HUMAN
Immunogen	
Immunogen	
Region	
Specificity	
Immunogen	
Sequence	