

Anti-CD366/HAVCR2 antibody [RMT3-23] {APC} (STJA0002099)

STJA0002099

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rat monoclonal antibody anti-Hepatitis A virus cellular receptor 2 is suitable for use in Flow Cytometry research applications.
Applications	FC
Host/Source	Rat
Reactivity	Mouse

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	RMT3-23
Concentration	0.2 mg/mL
Conjugation	APC
Purification	
Dilution Range	
Formulation	PBS with 0.05% Proclin300, 1% BSA
Isotype	IgG2ak
Storage Instruction	Recommend storing between 2-8°C and protecting from prolonged exposure to light. Do not freeze this product.

TARGET INFORMATION

Gene ID	84868
Gene Symbol	HAVCR2
Uniprot ID	HAVR2_HUMAN
Immunogen	
Immunogen Region	
Specificity	
Immunogen Sequence	