

Anti-CD25/IL2RA antibody [BC96] {PE/Cyanine5} (STJA0002108)

STJA0002108

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-Interleukin-2 receptor subunit alpha is suitable for use in Flow Cytometry research applications.
Applications	FC
Host/Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	BC96
Concentration	
Conjugation	PE/Cyanine5
Purification	
Dilution Range	
Formulation	PBS with 0.05% Proclin300, 1% BSA
Isotype	IgG1k
Storage Instruction	Recommend storing between 2-8°C and protecting from prolonged exposure to light. Do not freeze this product.

TARGET INFORMATION

Gene ID	3559
Gene Symbol	IL2RA
Uniprot ID	IL2RA_HUMAN
Immunogen	
Immunogen Region	
Specificity	
Immunogen Sequence	