

Anti-CD28 antibody [CD28.2] {AF488} (STJA0002118)

STJA0002118

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-T-cell-specific surface glycoprotein CD28 is suitable for use in Flow Cytometry research applications.
Applications	FC
Host/Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	CD28.2
Concentration	
Conjugation	AF488
Purification	
Dilution Range	
Formulation	PBS with 0.05% Proclin300, 1% BSA
Isotype	IgG1k
Storage Instruction	Recommend storing between 2-8°C and protecting from prolonged exposure to light. Do not freeze this product.

TARGET INFORMATION

Gene ID	940
Gene Symbol	CD28
Uniprot ID	CD28_HUMAN
Immunogen	
Immunogen Region	
Specificity	
Immunogen Sequence	