

## Anti-IFNG antibody [B27] (STJA0002121)

STJA0002121

### GENERAL INFORMATION

<b>Product Type</b>	Primary antibodies
<b>Short Description</b>	Mouse monoclonal antibody anti-Interferon gamma is suitable for use in IConfocal MicroscopyC, Western Blot, Immunohistochemistry and Immunoprecipitation research applications.
<b>Applications</b>	ICFC, WB, IHC, IP
<b>Host/Source</b>	Mouse
<b>Reactivity</b>	Human

### PRODUCT PROPERTIES

<b>Clonality</b>	Monoclonal
<b>Clone ID</b>	B27
<b>Concentration</b>	0.5 mg/mL
<b>Conjugation</b>	Unconjugated
<b>Purification</b>	
<b>Dilution Range</b>	
<b>Formulation</b>	PBS with 0.05% Proclin300
<b>Isotype</b>	IgG1k
<b>Storage</b>	Recommend storing at 2-8°C. Do not freeze this product.
<b>Instruction</b>	

### TARGET INFORMATION

<b>Gene ID</b>	<a href="#">3458</a>
<b>Gene Symbol</b>	<a href="#">IFNG</a>
<b>Uniprot ID</b>	<a href="#">IFNG_HUMAN</a>
<b>Immunogen</b>	
<b>Immunogen</b>	
<b>Region</b>	
<b>Specificity</b>	
<b>Immunogen</b>	
<b>Sequence</b>	