

Anti-IL2 antibody [JES6-5H4] {APC} (STJA0002146)

STJA0002146

GENERAL INFORMATION


| | |
|--------------------------|--|
| Product Type | Primary antibodies |
| Short Description | Rat monoclonal antibody anti-Interleukin-2 is suitable for use in IConfocal MicroscopyC research applications. |
| Applications | ICFC |
| Host/Source | Rat |
| Reactivity | Mouse |

PRODUCT PROPERTIES

| | |
|----------------------------|--|
| Clonality | Monoclonal |
| Clone ID | JES6-5H4 |
| Concentration | |
| Conjugation | APC |
| Purification | |
| Dilution Range | |
| Formulation | PBS with 0.05% Proclin300, 1% BSA |
| Isotype | IgG2bk |
| Storage Instruction | Recommend storing between 2-8°C and protecting from prolonged exposure to light. Do not freeze this product. |

TARGET INFORMATION

| | |
|---------------------------|---------------------------|
| Gene ID | 16183 |
| Gene Symbol | IL2 |
| Uniprot ID | IL2_MOUSE |
| Immunogen | |
| Immunogen Region | |
| Specificity | |
| Immunogen Sequence | |



This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081