

## Anti-IL4 antibody [MP4-25D2] (STJA0002153)

STJA0002153

### GENERAL INFORMATION

<b>Product Type</b>	Primary antibodies
<b>Short Description</b>	Rat monoclonal antibody anti-Interleukin-4 is suitable for use in IConfocal MicroscopyC, ELISA and Immunohistochemistry research applications.
<b>Applications</b>	ICFC, ELISA, IHC
<b>Host/Source</b>	Rat
<b>Reactivity</b>	Human

### PRODUCT PROPERTIES

<b>Clonality</b>	Monoclonal
<b>Clone ID</b>	MP4-25D2
<b>Concentration</b>	0.5 mg/mL
<b>Conjugation</b>	Unconjugated
<b>Purification</b>	
<b>Dilution Range</b>	
<b>Formulation</b>	PBS with 0.05% Proclin300
<b>Isotype</b>	IgG1k
<b>Storage Instruction</b>	Recommend storing at 2-8°C. Do not freeze this product.

### TARGET INFORMATION

<b>Gene ID</b>	<a href="#">3565</a>
<b>Gene Symbol</b>	<a href="#">IL4</a>
<b>Uniprot ID</b>	<a href="#">IL4_HUMAN</a>
<b>Immunogen</b>	
<b>Immunogen</b>	
<b>Region</b>	
<b>Specificity</b>	
<b>Immunogen</b>	
<b>Sequence</b>	