

## Anti-IL4 antibody [MP4-25D2] {PE} (STJA0002155)

STJA0002155

### GENERAL INFORMATION

<b>Product Type</b>	Primary antibodies
<b>Short Description</b>	Rat monoclonal antibody anti-Interleukin-4 is suitable for use in IConfocal MicroscopyC research applications.
<b>Applications</b>	ICFC
<b>Host/Source</b>	Rat
<b>Reactivity</b>	Human

### PRODUCT PROPERTIES

<b>Clonality</b>	Monoclonal
<b>Clone ID</b>	MP4-25D2
<b>Concentration</b>	
<b>Conjugation</b>	PE
<b>Purification</b>	
<b>Dilution Range</b>	
<b>Formulation</b>	PBS with 0.05% Proclin300, 1% BSA
<b>Isotype</b>	IgG1k
<b>Storage Instruction</b>	Recommend storing between 2-8°C and protecting from prolonged exposure to light. Do not freeze this product.

### TARGET INFORMATION

<b>Gene ID</b>	<a href="#">3565</a>
<b>Gene Symbol</b>	<a href="#">IL4</a>
<b>Uniprot ID</b>	<a href="#">IL4_HUMAN</a>
<b>Immunogen</b>	
<b>Immunogen Region</b>	
<b>Specificity</b>	
<b>Immunogen Sequence</b>	