

Anti-IL4 antibody [11B11] (STJA0002157)

STJA0002157

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rat monoclonal antibody anti-Interleukin-4 is suitable for use in Universal, Universal, IConfocal MicroscopyC, Western Blot, Immunohistochemistry and Immunoprecipitation research applications.
Applications	Neut, Stim, ICFC, WB, IHC, IP
Host/Source	Rat
Reactivity	Mouse

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	11B11
Concentration	0.5 mg/mL
Conjugation	Unconjugated
Purification	
Dilution Range	
Formulation	0.2 µm filtered in PBS, pH7.2. Azide Free (AF) /Low Endotoxin (LE) : Contains no stabilizers or preservatives. Endotoxin level is < 2 EU/mg as determined by LAL gel clotting assay.
Isotype	IgG1k
Storage	Recommend storing at 2-8°C. Do not freeze this product.
Instruction	

TARGET INFORMATION

Gene ID	16189
Gene Symbol	IL4
Uniprot ID	IL4_MOUSE
Immunogen	
Immunogen Region	
Specificity	
Immunogen Sequence	