

Anti-IL4 antibody [11B11] (STJA0002158)

STJA0002158

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rat monoclonal antibody anti-Interleukin-4 is suitable for use in IConfocal MicroscopyC, ELISA, Immunohistochemistry and Immunoprecipitation research applications.
Applications	ICFC, ELISA, IHC, IP
Host/Source	Rat
Reactivity	Mouse

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	11B11
Concentration	0.5 mg/mL
Conjugation	Unconjugated
Purification	
Dilution Range	
Formulation	PBS with 0.05% Proclin300
Isotype	IgG1k
Storage	Recommend storing at 2-8°C. Do not freeze this product.
Instruction	

TARGET INFORMATION

Gene ID	16189
Gene Symbol	IL4
Uniprot ID	IL4_MOUSE
Immunogen	
Immunogen	
Region	
Specificity	
Immunogen	
Sequence	