

Anti-IL-5 antibody [TRFK5] (STJA0002163)

STJA0002163

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rat monoclonal antibody anti-IL-5 is suitable for use in Universal, Universal, IConfocal MicroscopyC, Western Blot and Immunohistochemistry research applications.
Applications	Depletion, Neut, ICFC, WB, IHC
Host/Source	Rat
Reactivity	Human, Mouse

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	TRFK5
Concentration	0.5 mg/mL
Conjugation	Unconjugated
Purification	
Dilution Range	
Formulation	0.2 µm filtered in PBS, pH7.2. Azide Free (AF) /Low Endotoxin (LE) : Contains no stabilizers or preservatives. Endotoxin level is < 2 EU/mg as determined by LAL gel clotting assay.
Isotype	IgG1k
Storage Instruction	Recommend storing at 2-8°C. Do not freeze this product.

TARGET INFORMATION

Gene ID	3567
Gene Symbol	IL5
Uniprot ID	IL5_HUMAN
Immunogen	
Immunogen	
Region	
Specificity	
Immunogen	
Sequence	