

## Anti-IL6 antibody [MP5-20F3] (STJA0002173)

STJA0002173

### GENERAL INFORMATION

|                          |  |
|--------------------------|--|
| <b>Product Type</b>      | Primary antibodies   |
| <b>Short Description</b> | Rat monoclonal antibody anti-Interleukin-6 is suitable for use in IConfocal MicroscopyC, ELISA and Immunohistochemistry research applications. |
| <b>Applications</b>      | ICFC, ELISA, IHC   |
| <b>Host/Source</b>       | Rat  |
| <b>Reactivity</b>        | Mouse  |

### PRODUCT PROPERTIES

|                            |   |
|----------------------------|---|
| <b>Clonality</b>           | Monoclonal  |
| <b>Clone ID</b>            | MP5-20F3  |
| <b>Concentration</b>       | 0.5 mg/mL   |
| <b>Conjugation</b>         | Unconjugated  |
| <b>Purification</b>        |   |
| <b>Dilution Range</b>      |   |
| <b>Formulation</b>         | PBS with 0.05% Proclin300                               |
| <b>Isotype</b>             | IgG1k   |
| <b>Storage Instruction</b> | Recommend storing at 2-8°C. Do not freeze this product. |

### TARGET INFORMATION

|                    |                           |
|--------------------|---------------------------|
| <b>Gene ID</b>     | <a href="#">16193</a>     |
| <b>Gene Symbol</b> | <a href="#">IL6</a>       |
| <b>Uniprot ID</b>  | <a href="#">IL6_MOUSE</a> |
| <b>Immunogen</b>   |                           |
| <b>Immunogen</b>   |                           |
| <b>Region</b>      |                           |
| <b>Specificity</b> |                           |
| <b>Immunogen</b>   |                           |
| <b>Sequence</b>    |                           |