

## Anti-CD20/MS4A1 antibody [2H7] (STJA0002192)

STJA0002192

### GENERAL INFORMATION

<b>Product Type</b>	Primary antibodies
<b>Short Description</b>	Mouse monoclonal antibody anti-B-lymphocyte antigen CD20 is suitable for use in Universal, Flow Cytometry, Immunohistochemistry and Immunoprecipitation research applications.
<b>Applications</b>	Depletion, FC, IHC, IP
<b>Host/Source</b>	Mouse
<b>Reactivity</b>	Human

### PRODUCT PROPERTIES

<b>Clonality</b>	Monoclonal
<b>Clone ID</b>	2H7
<b>Concentration</b>	0.5 mg/mL
<b>Conjugation</b>	Unconjugated
<b>Purification</b>	
<b>Dilution Range</b>	
<b>Formulation</b>	0.2 µm filtered in PBS, pH7.2. Azide Free (AF) /Low Endotoxin (LE) : Contains no stabilizers or preservatives. Endotoxin level is < 2 EU/mg as determined by LAL gel clotting assay.
<b>Isotype</b>	IgG2bk
<b>Storage</b>	Recommend storing at 2-8°C. Do not freeze this product.
<b>Instruction</b>	

### TARGET INFORMATION

<b>Gene ID</b>	931
<b>Gene Symbol</b>	MS4A1
<b>Uniprot ID</b>	CD20_HUMAN
<b>Immunogen</b>	
<b>Immunogen</b>	
<b>Region</b>	
<b>Specificity</b>	
<b>Immunogen</b>	
<b>Sequence</b>	