

Anti-PDCD1 antibody [J110] (STJA0002201)

STJA0002201

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-Programmed cell death protein 1 is suitable for use in Flow Cytometry and Immunohistochemistry research applications.
Applications	FC, IHC
Host/Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	J110
Concentration	0.5 mg/mL
Conjugation	Unconjugated
Purification	
Dilution Range	
Formulation	PBS with 0.05% Proclin300
Isotype	IgG1k
Storage	Recommend storing at 2-8°C. Do not freeze this product.
Instruction	

TARGET INFORMATION

Gene ID	5133
Gene Symbol	PDCD1
Uniprot ID	PDCD1_HUMAN
Immunogen	
Immunogen	
Region	
Specificity	
Immunogen	
Sequence	