

Anti-CD40 antibody [G28.5] {PE/Cyanine5} (STJA0002204)

STJA0002204

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-Tumor necrosis factor receptor superfamily member 5 is suitable for use in Flow Cytometry research applications.
Applications	FC
Host/Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	G28.5
Concentration	
Conjugation	PE/Cyanine5
Purification	
Dilution Range	
Formulation	PBS with 0.05% Proclin300, 1% BSA
Isotype	IgG1k
Storage	Recommend storing between 2-8°C and protecting from prolonged exposure to light. Do not freeze this product.
Instruction	

TARGET INFORMATION

Gene ID	958
Gene Symbol	CD40
Uniprot ID	TNR5_HUMAN
Immunogen	
Immunogen	
Region	
Specificity	
Immunogen	
Sequence	