

Anti-H-2 antibody [M1/42] (STJA0002215)

STJA0002215

GENERAL INFORMATION

| | |
|--------------------------|---|
| Product Type | Primary antibodies |
| Short Description | Rat monoclonal antibody anti-H-2 is suitable for use in Flow Cytometry and Immunoprecipitation research applications. |
| Applications | FC, IP |
| Host/Source | Rat |
| Reactivity | Mouse |

PRODUCT PROPERTIES

| | |
|----------------------------|---|
| Clonality | Monoclonal |
| Clone ID | M1/42 |
| Concentration | 0.5 mg/mL |
| Conjugation | Unconjugated |
| Purification | |
| Dilution Range | |
| Formulation | PBS with 0.05% Proclin300 |
| Isotype | IgG2ak |
| Storage Instruction | Recommend storing at 2-8°C. Do not freeze this product. |

TARGET INFORMATION

| | |
|---------------------------|--|
| Gene ID | |
| Gene Symbol | |
| Uniprot ID | |
| Immunogen | |
| Immunogen Region | |
| Specificity | |
| Immunogen Sequence | |