

Anti-CD16/CD80 antibody [3G8] {PerCP/Cyanine5.5} (STJA0002248)

STJA0002248

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-T-lymphocyte activation antigen CD80 is suitable for use in Flow Cytometry research applications.
Applications	FC
Host/Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	3G8
Concentration	
Conjugation	PerCP/Cyanine5.5
Purification	
Dilution Range	
Formulation	PBS with 0.05% Proclin300, 1% BSA
Isotype	IgG1k
Storage Instruction	Recommend storing between 2-8°C and protecting from prolonged exposure to light. Do not freeze this product.

TARGET INFORMATION

Gene ID	941
Gene Symbol	CD80
Uniprot ID	CD80_HUMAN
Immunogen	
Immunogen Region	
Specificity	
Immunogen	
Sequence	