

Anti-IgG2ak antibody [C1.18.4] (STJA0002370)

STJA0002370

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-IgG2a, k is suitable for use in Flow Cytometry research applications.
Applications	FC
Host/Source	Mouse
Reactivity	

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	C1.18.4
Concentration	0.5 mg/mL
Conjugation	Unconjugated
Purification	
Dilution Range	
Formulation	0.2 µm filtered in PBS, pH7.2. Azide Free (AF) /Low Endotoxin (LE) : Contains no stabilizers or preservatives. Endotoxin level is < 2 EU/mg as determined by LAL gel clotting assay.
Isotype	IgG2ak
Storage Instruction	Recommend storing at 2-8°C. Do not freeze this product.

TARGET INFORMATION

Gene ID
Gene Symbol
Uniprot ID
Immunogen
Immunogen Region
Specificity
Immunogen Sequence