

Anti-IgG2ak antibody [2A3] {PE/Cyanine5} (STJA0002438)

STJA0002438

GENERAL INFORMATION

| | |
|--------------------------|----------------------------------------------------------------------------------------------------|
| Product Type | Primary antibodies |
| Short Description | Rat monoclonal antibody anti-IgG2a, k is suitable for use in Flow Cytometry research applications. |
| Applications | FC |
| Host/Source | Rat |
| Reactivity | |

PRODUCT PROPERTIES

| | |
|----------------------------|--------------------------------------------------------------------------------------------------------------|
| Clonality | Monoclonal |
| Clone ID | 2A3 |
| Concentration | |
| Conjugation | PE/Cyanine5 |
| Purification | |
| Dilution Range | |
| Formulation | PBS with 0.05% Proclin300, 1% BSA |
| Isotype | IgG2ak |
| Storage Instruction | Recommend storing between 2-8°C and protecting from prolonged exposure to light. Do not freeze this product. |

TARGET INFORMATION

| | |
|---------------------------|--|
| Gene ID | |
| Gene Symbol | |
| Uniprot ID | |
| Immunogen | |
| Immunogen Region | |
| Specificity | |
| Immunogen Sequence | |