

Anti-IgG antibody [PIP] {PerCP/Cyanine5.5} (STJA0002491)

STJA0002491

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Armenian Hamster monoclonal antibody anti-IgG is suitable for use in Flow Cytometry research applications.
Applications	FC
Host/Source	Armenian Hamster
Reactivity	

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	PIP
Concentration	
Conjugation	PerCP/Cyanine5.5
Purification	
Dilution Range	
Formulation	PBS with 0.05% Proclin300, 1% BSA
Isotype	IgG
Storage Instruction	Recommend storing between 2-8°C and protecting from prolonged exposure to light. Do not freeze this product.

TARGET INFORMATION

Gene ID
Gene Symbol
Uniprot ID
Immunogen
Immunogen Region
Specificity
Immunogen Sequence