

Anti-CD3 antibody [OKT3, 3G8, 5.1H11] {FITC/PE} (STJA0002509)

STJA0002509

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	NA monoclonal antibody anti-CD3-FITC/CD16+CD56-PE is suitable for use in Flow Cytometry research applications.
Applications	FC
Host/Source	
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	OKT3, 3G8, 5.1H11
Concentration	
Conjugation	FITC/PE
Purification	
Dilution Range	
Formulation	PBS with 0.05% Proclin300, 1% BSA
Isotype	
Storage Instruction	Recommend storing between 2-8°C and protecting from prolonged exposure to light. Do not freeze this product.

TARGET INFORMATION

Gene ID	
Gene Symbol	
Uniprot ID	
Immunogen	
Immunogen Region	
Specificity	
Immunogen Sequence	