

Anti-CD3 antibody [OKT3, RPA-T4, HI30] {FITC/PE/PerCP-Cyanine5.5} (STJA0002511)

STJA0002511

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	NA monoclonal antibody anti-CD3-FITC/CD4-PE/CD45-PerCP-Cyanine5.5 is suitable for use in Flow Cytometry research applications.
Applications	FC
Host/Source	
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	OKT3, RPA-T4, HI30
Concentration	
Conjugation	FITC/PE/PerCP-Cyanine5.5
Purification	
Dilution Range	
Formulation	PBS with 0.05% Proclin300, 1% BSA
Isotype	
Storage Instruction	Recommend storing between 2-8°C and protecting from prolonged exposure to light. Do not freeze this product.

TARGET INFORMATION

Gene ID	
Gene Symbol	
Uniprot ID	
Immunogen	
Immunogen Region	
Specificity	
Immunogen Sequence	