

**Anti-IgG2a, k-antibody [C1.18.4] {PE/TR} (STJA0002685)**  
STJA0002685**GENERAL INFORMATION**

<b>Product Type</b>	Primary antibodies
<b>Short Description</b>	Mouse monoclonal antibody anti-IgG2a, k is suitable for use in Flow Cytometry.
<b>Applications</b>	FC
<b>Host/Source</b>	Mouse
<b>Reactivity</b>	

**PRODUCT PROPERTIES**

<b>Clonality</b>	Monoclonal
<b>Clone ID</b>	C1.18.4
<b>Concentration</b>	
<b>Conjugation</b>	PE/TR
<b>Purification</b>	
<b>Dilution Range</b>	
<b>Formulation</b>	PBS with 0.05% Proclin300, 1% BSA
<b>Isotype</b>	IgG2a, k
<b>Storage Instruction</b>	Store at 2~8°C and protected from prolonged exposure to light. Do not freeze.

**TARGET INFORMATION**

<b>Gene ID</b>	
<b>Gene Symbol</b>	
<b>Uniprot ID</b>	
<b>Immunogen</b>	
<b>Immunogen Region</b>	
<b>Specificity</b>	
<b>Immunogen Sequence</b>	