

Anti-IgG2b, k-antibody [MPC-11] {PE/TR} (STJA0002687)

STJA0002687

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-IgG2b, k is suitable for use in Flow Cytometry.
Applications	FC
Host/Source	Mouse
Reactivity	

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	MPC-11
Concentration	
Conjugation	PE/TR
Purification	
Dilution Range	
Formulation	PBS with 0.05% Proclin300, 1% BSA
Isotype	IgG2b, k
Storage Instruction	Store at 2~8°C and protected from prolonged exposure to light. Do not freeze.

TARGET INFORMATION

Gene ID	
Gene Symbol	
Uniprot ID	
Immunogen	
Immunogen Region	
Specificity	
Immunogen Sequence	