

Anti-IgG-antibody [PIP] {PE/TR} (STJA0002695)

STJA0002695

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Armenian Hamster monoclonal antibody anti-IgG is suitable for use in Flow Cytometry.
Applications	FC
Host/Source	Armenian Hamster
Reactivity	

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	PIP
Concentration	
Conjugation	PE/TR
Purification	
Dilution Range	
Formulation	PBS with 0.05% Proclin300, 1% BSA
Isotype	IgG
Storage Instruction	Store at 2~8°C and protected from prolonged exposure to light. Do not freeze.

TARGET INFORMATION

Gene ID	
Gene Symbol	
Uniprot ID	
Immunogen	
Immunogen Region	
Specificity	
Immunogen Sequence	