

**Anti-IgG2a, Lambda antibody [B39-4] {PE/Cyanine5.5} (STJA0002757)**

STJA0002757

**GENERAL INFORMATION**

|                          |   |
|--------------------------|---|
| <b>Product Type</b>      | Primary antibodies  |
| <b>Short Description</b> | Rat monoclonal antibody anti-IgG2a, Lambda is suitable for use in Flow Cytometry research applications. |
| <b>Applications</b>      | FC  |
| <b>Host/Source</b>       | Rat   |
| <b>Reactivity</b>        |   |

**PRODUCT PROPERTIES**

|                            |  |
|----------------------------|--|
| <b>Clonality</b>           | Monoclonal   |
| <b>Clone ID</b>            | B39-4  |
| <b>Concentration</b>       |  |
| <b>Conjugation</b>         | PE/Cyanine5.5  |
| <b>Purification</b>        |  |
| <b>Dilution Range</b>      |  |
| <b>Formulation</b>         | PBS with 0.05% Proclin300, 1% BSA  |
| <b>Isotype</b>             | IgG2a, lambda  |
| <b>Storage Instruction</b> | Recommend storing between 2-8°C and protecting from prolonged exposure to light. Do not freeze this product. |

**TARGET INFORMATION**

|                           |  |
|---------------------------|--|
| <b>Gene ID</b>            |  |
| <b>Gene Symbol</b>        |  |
| <b>Uniprot ID</b>         |  |
| <b>Immunogen</b>          |  |
| <b>Immunogen Region</b>   |  |
| <b>Specificity</b>        |  |
| <b>Immunogen Sequence</b> |  |