

Anti-IgG2a, Lambda-antibody [B39-4] {PE/TR} (STJA0002762)
STJA0002762**GENERAL INFORMATION**

Product Type	Primary antibodies
Short Description	Rat monoclonal antibody anti-IgG2a, Lambda is suitable for use in Flow Cytometry.
Applications	FC
Host/Source	Rat
Reactivity	

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	B39-4
Concentration	
Conjugation	PE/TR
Purification	
Dilution Range	
Formulation	PBS with 0.05% Proclin300, 1% BSA
Isotype	IgG2a, lambda
Storage Instruction	Store at 2~8°C and protected from prolonged exposure to light. Do not freeze.

TARGET INFORMATION

Gene ID	
Gene Symbol	
Uniprot ID	
Immunogen	
Immunogen Region	
Specificity	
Immunogen Sequence	