

Anti-IgG2a, Lambda antibody [B39-4] (STJA0002764)

STJA0002764

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rat monoclonal antibody anti-IgG2a, Lambda is suitable for use in Flow Cytometry research applications.
Applications	FC
Host/Source	Rat
Reactivity	

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	B39-4
Concentration	0.5 mg/mL
Conjugation	Unconjugated
Purification	
Dilution Range	
Formulation	0.2 µm filtered in PBS, pH7.2. Contains no stabilizers or preservatives. Endotoxin level is < 2 EU/mg (<0.002EU/µg) as determined by LAL gel clotting assay. Purity is >95% as determined by SDS-PAGE.
Isotype	IgG2a, lambda
Storage Instruction	Recommend storing at 2-8°C. Do not freeze this product.

TARGET INFORMATION

Gene ID	
Gene Symbol	
Uniprot ID	
Immunogen	
Immunogen Region	
Specificity	
Immunogen Sequence	