

Anti-DNase I antibody (STJA0002809)

STJA0002809

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-DNase I is suitable for use in ELISA and Western Blot research applications.
Applications	ELISA, WB
Host/Source	Rabbit
Reactivity	DNASE1

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	
Conjugation	Unconjugated
Purification	Bovine DNaseI-affinity purification
Dilution Range	ELISA 0.05 to 0.1 µg/ml; WB 0.1 to 0.2 µg/ml; IF 0.2 to 12 Mu g/ml Antibody working titer has to be established practically for each particular antigen and assay format.
Formulation	Concentrated ammonium sulphate in PBS pH 7.4.
Isotype	IgG
Storage Instruction	Store at +4°C upon receipt. As product is (NH ₄) ₂ SO ₄ precipitate, mix well by pipetting or vortexing prior use. Do not freeze

TARGET INFORMATION

Gene ID
Gene Symbol
Uniprot ID
Immunogen
Immunogen
Region
Specificity
Immunogen
Sequence