

Anti-HIV-1 Tat antibody (59-67aa) [N8] (STJA0002850)

STJA0002850

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-HIV-1 Tat (59-67aa) is suitable for use in ELISA, Western Blot and ChImmunoprecipitation research applications.
Applications	ELISA, WB, ChIP
Host/Source	Mouse
Reactivity	HIV-1

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	N8
Concentration	
Conjugation	Unconjugated
Purification	Protein G purification
Dilution Range	ELISA 0.05 to 0.1 µg/ml; WB 2 µg/ml. Monoclonal antibody working titer has to be established practically for each particular antigen and assay format
Formulation	PBS pH 7.4, with 0.1% sodium azide
Isotype	IgG2b
Storage	Store at -20 to -70 °C upon receipt. Divide antibody into aliquots prior usage. Avoid multiple freeze-thaw cycles
Instruction	

TARGET INFORMATION

Gene ID	
Gene Symbol	
Uniprot ID	
Immunogen	Recombinant Tat protein of HIV-1 HAN (purified from E. coli). Mapped to amino acids 59-67 (PQDSQTHQV – B subtype a.a. consensus)
Immunogen Region	59-67aa
Specificity	Reacts in ELISA and WB with HIV-1 Tat of HXB2, BH10, IIB and HAN strains. Mapped to amino acids 59-67 (PQDSQTHQV – B subtype a.a. consensus)
Immunogen Sequence	

