

Anti-MLV p30 antibody (STJA0002864)

STJA0002864

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-MLV p30 is suitable for use in ELISA and Western Blot research applications.
Applications	ELISA, WB
Host/Source	Rabbit
Reactivity	MLV

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	
Conjugation	Unconjugated
Purification	MLV p30 affinity purified
Dilution Range	ELISA 0.01 ug/ml; WB 0.5-1 µg/ml Polyclonal antibody working titer has to be established practically for each particular antigen and assay format
Formulation	Concentrated ammonium sulphate in PBS pH 7.4.
Isotype	IgG
Storage Instruction	Store at +4°C upon receipt. As product is (NH ₄) ₂ SO ₄ precipitate, mix well by pipetting or vortexing prior use

TARGET INFORMATION

Gene ID	
Gene Symbol	
Uniprot ID	
Immunogen	Full length recombinant MLV-Gag p30 capsid protein, isolate FB29, purified from E. coli. MLV p30
Immunogen Region	
Specificity	MLV Gag, p30
Immunogen Sequence	