

## Anti-Hev b3 antibody [5] (STJA0002870)

STJA0002870

### GENERAL INFORMATION

<b>Product Type</b>	Primary antibodies
<b>Short Description</b>	Mouse monoclonal antibody anti-Small rubber particle protein is suitable for use in ELISA and Western Blot research applications.
<b>Applications</b>	ELISA, WB
<b>Host/Source</b>	Mouse
<b>Reactivity</b>	HEVB3

### PRODUCT PROPERTIES

<b>Clonality</b>	Monoclonal
<b>Clone ID</b>	5
<b>Concentration</b>	
<b>Conjugation</b>	Unconjugated
<b>Purification</b>	Protein G purification
<b>Dilution Range</b>	ELISA 14 000. Reacts as binding antibody in capture ELISA with Hev b3 detection monoclonal antibody A3-705-100. WB 11000 Monoclonal antibody working titer has to be established practically for each particular antigen and assay format
<b>Formulation</b>	PBS, 0.1% sodium azide
<b>Isotype</b>	IgG1
<b>Storage Instruction</b>	Store at -20 to -70 °C upon receipt. Divide antibody into aliquots prior usage. Avoid multiple freeze-thaw cycles as product degradation may result

### TARGET INFORMATION

<b>Gene ID</b>	
<b>Gene Symbol</b>	
<b>Uniprot ID</b>	
<b>Immunogen Region</b>	Hev b3-MBP fusion protein expressed and purified from E. coli. Hev b3
<b>Specificity</b>	Recombinant and native Hev b3
<b>Immunogen Sequence</b>	